



PRODUCT CATALOGUE

2024/2025

Dear Esteemed Customers & Partners

We express our sincere thanks and appreciate your business association with us since last so many years. Despite a lot of choices and alternatives available for you, you kept your faith on us. This trust on us keeps us motivated to serve you all more passionately. Whenever you are looking for something special and trustworthy in IVD, we hope that we shall be your first choice in this advanced diagnostic field.

Our goal is not only to be a leader in diagnostic industry by developing wide range of Biochemistry, Hematology, Rapid, Elisa & Molecular kits but also to excel in taking care of human health more closely and carefully through POC products.

Our goal remains to be focussed on introduction of innovative and high-tech products for providing better diagnostic service through make in India initiative. With emphasis on Q-Line Biotech playing a vital role towards prevention & diagnosis of diseases and making it a responsible and reliable global supplier of IVD products, we once again thank our all associates.

Thank you all.

About Company

Q-Line Biotech is an IVD company based in national capital of India, New Delhi. Earlier known as POCT Services Pvt. Ltd. company started its journey in 2014 taking over the business of Merck Diagnostics.

Presently company is involved in manufacturing and distribution of various Clinical Chemistry Reagents, Rapid Tests, Elisa Kits, POC products and Molecular kits as its core business. Company focusses on hard core innovations in Diagnostics, Molecular, Bio processing and various IVD analyzers. Q-Line Biotech has joint ventures with various international companies for manufacturing of POC products and IVD analyzers for lab use.

Establishing highest quality standards to meet the growing market demands and challenges, we are an ISO 9001:2015 certified company with well established R & D facility and state of the art GMP & ISO :13486:2016 certified manufacturing facility at New Delhi.

Presently we are one of the fastest growing and expanding Indian IVD company with market presence in Indian subcontinent and overseas.

With proud owner of "Q-Line" brand our vision is to transpire Q-LINE BIOTECH as an eminent player in global market for all healthcare requirements & providing a Complete Hospital Solution to the Medical fraternity. With complete focus on excellence in terms of standards of quality & affordability.



TABLE OF CONTENTS

Semi Automatic Biochemistry Analyzers	1-1	System Pack Reagents	19-21
Fully Automatic Biochemistry Analyzers	2-3	Special Chemistry Reagents	21-22
Hematology Analyzers	4-5	Diazyme Reagents	22-22
ESR Analyzers	5-5	Coagulation Reagents	23-23
Hemoglobin Testing	5-5	Electrolyte Reagents	23-23
Electrolyte Analyzers	5-6	Urine Chemistry Strips	24-24
Urine Chemistry Analyzers	7-8	Rapid Test Range	24-24
Elisa Readers & Washers	8-10	Elisa Kits	24-24
Biomedical Systems	10-12	Molecular Range	25-25
Fully Automated Nucleic Acid Extraction Systems	13-13	Hematology Reagents, Calibrator, Controls & ESR Tubes	25-25
Reagents & Consumables	14-26	P.O.C. Products	26-26
General Biochemistry Reagents	15-19	Elitech Biomedical Reagents	26-26

Semi Automatic Biochemistry Analyzers

Microlab 300



Microlab 300

- Maintenance free bellows pump technology.
- Superior optical design gives accuracy and better filter life.
- High quality inbuilt QWERTY keyboard & external keyboard facility.
- Multiple replicates measurement facility of samples & calibrator for better accuracy.
- User programmable non-linearity % for hyperactive sample verification.
- Seven standard filters.
- Excellent performance & rugged hardware.
- USFDA Listed & CE Certified.

Ordering Information:

Product Code	Product Name
6003075VS	Microlab 300 LX (w/o inbuilt printer)
6002310VS	Microlab 300 EX (with inbuilt printer)

Innolab 200

- Windows based user friendly interface.
- 7.0 Inch touch screen LCD display.
- Online reaction curve.
- Up to 300 tests & 2,00,000 patient results storage.
- Auto sleep function for lamp protection & cuvette mode.
- In built thermal printer, external keyboard and mouse attachable.
- 5 USB ports with ethernet port for LIS.
- USFDA Listed & CE certified




InnoLab
200

Ordering Information:

Product Code	Product Name
600002PS	Innolab 200

Fully Automatic Biochemistry Analyzers

Selectra Pro M Lite



SELECTRA PRO M LITE
CHEMISTRY SYSTEM

- Random access clinical chemistry analyzer.
- Throughput - Up to 180 test /hr.
- On board laundry and reagent cooling.
- Separate rotors with 32 reagent and 50 sample positions.
- High precision separate syringes for sample & reagents.
- Separate reagent & sample probes with level detection.
- Integrated mixer and collision protection.
- Least water consumption of ~ 500 mL /hr.
- Cost effective, minimum manpower.
- USFDA Listed & CE Certified

Ordering Information:

Product Code	Product Name
6003439VS	Selectra Pro M LITE

Selectra Pro S Lite

- Random access clinical chemistry analyser.
- Throughput - Up to 133 test /hr.
- On board laundry and reagent cooling.
- Combined reagent/sample rotor with 30 reagent and 25 sample positions.
- High precision separate syringes for sample & reagents.
- Reagent / sample probe with level detector.
- Integrated mixer and collision protection.
- Water consumption of ~ 950 mL /hr.
- USFDA Listed & CE certified



SELECTRA PRO S LITE
CHEMISTRY SYSTEM

Ordering Information:

Product Code	Product Name
6003548VS	Selectra Pro S LITE

Fully Automatic Biochemistry Analyzers

ASTIRA C-MS



- Random access clinical chemistry analyzer.
- Throughput constant 240 tests/hr.
- On board laundry and reagent cooling.
- Grating photometry technology.
- Multifunctional carousel with 80 reagent and 40 sample positions flexible.
- In built sample & reagent barcode reader.
- Wavelength measuring range from 340-800 nm (12 wavelengths).
- Barcoded reagent packs.
- Long life high precision ceramic piston syringes.
- USFDA Listed & CE certified

Ordering Information:

Product Code	Product Name
ACMS0240QB	Astira C-MS

ASTIRA CX400

- Random access clinical chemistry analyzer.
- Throughput 400 T/H (constant).
- On board laundry and reagent cooling.
- Grating photometry technology.
- 115 sample positions, 112 reagent positions.
- 7 stops 11 steps rising systems.
- In built sample & reagent barcode reader.
- Can handle up-to 4 reagent chemistry.
- Wavelength measuring range from 340-800 nm (12 wavelengths).
- Barcoded reagent packs.
- Liquid level detection, clot detection & collision protection.
- LIS interface, support the bi-directional communication.
- USFDA Listed & CE certified.



Ordering Information:

Product Code	Product Name
ACCX0400QB	Astira CX400

Hemoglobin Analyzer

HbA1c (HPLC)



Q-A1c H10

- Automated system based on High Pressure Liquid Chromatography (HPLC).
- High Precision with reliable column, $CV \leq 1.5\%$.
- Variant recognition with β -Thalassemia mode.
- Auto sample station with 10 positions with barcode scanning (external).
- Primary tube sampling with cap piercing.
- Variant Mode – Reported parameters HbA1a, HbA1b, HbF, LA1c, HbA1c, HbA0, V-Window.
- Thalassemia Mode- Reported parameters HbA1a, HbA1b, HbF, LA1c, HbA1c, HbA0, HbA2, HbE/HbD/HbS/HbC.
- HbA1c results reported in IFCC & NGSP values with ADAG value.
- Dual wavelength detection to avoid interference.
- Inbuilt thermal printer & external laser printer.
- QC Curve.
- USB & LAN, Ports, LIS Compatible.

Ordering Information:

Product Code	Product Name
QA1CH10001QB	Q-A1c H10 Hemoglobin Analyzer

HbA1c (HPLC)

- Automated system based on High Pressure Liquid Chromatography (HPLC).
- HbA1c Results within 130 Seconds.
- Reports HbA1c, HbF & Variant Window.
- High Precision with reliable column, $CV \leq 1.5\%$.
- Auto sample station with 5 samples loading capacity, suitable for medium to small labs.
- Primary tube sampling with barcode scanning (external).
- 10.1" TFT True Colour LCD Touch Screen Linux Software.
- HbA1c results reported in IFCC & NGSP values with ADAG value.
- Dual wavelength detection to avoid interference.
- Inbuilt thermal printer & external laser printer option.
- QC Curve.
- USB & LAN Ports, LIS Compatible.

COMING
SOON

Q-A1c H5

Ordering Information:

Product Code	Product Name
QA1CH50001QB	Q-A1c H5 Hemoglobin Analyzer

Haematology Analyzers

Medonic M32 Series - 3 Part Diff



Medonic

- 22 parameter three part diff analyzer.
- Throughput - 60 samples /hr.
- Based on high precision SRV technology.
- Floating discriminator & extended PLT count technology.
- Two reagent system with no mandatory start-up & shut down cycles.
- MPA for paediatric & donor screening from finger prick.
- High linearity PLT counting mode for blood bank application.
- Total quality concept with own Cal & QC .
- In-built blood mixer for uniform sample mixing .
- 1 USB device/4 USB host/1 LAN port .
- CE certified

Ordering Information:

Product Code	Product Name
1420021BMR	Medonic M32B
1420022BMR	M32M (With MPA device)



Medonic

Medonic M51 - 5 Part Diff

- 29 parameter Five-part diff analyzer with AL & IG.
- Throughput - Up to 60 samples/ hr.
- Tri-angle laser scatter Flow Cytometry (FCM) method.
- Four 2D scatter grams and three histograms.
- Open tube sample aspiration.
- Three reagents & cleaner for 5 diff analysis.
- Total quality concept with own Cal & QC.
- Support USB, LAN for bi-directional LIS.
- 50,000 patient sample memory.
- CE certified.



Medonic

Ordering Information:

Product Code	Product Name
1620020BM	Medonic M51

ESR Analyzers

Qualcyte® 10



Qualcyte® 10

- Random access ESR analyser.
- Based on Infrared laser beam technology.
- 10 Channels.
- Throughput 20 samples /hr.
- Reading time: 30 minutes.
- Excellent correlation with the manual Westergren method.
- External thermal printer (optional).
- LIS Interface & sample bar code reader (optional).
- Up-to 10000 samples storage.
- CE certified

Ordering Information:

Product Code	Product Name
SPXE10SG	Qualcyte 10
SPXE10PSG	Qualcyte10plus

Hemoglobin Testing System

Smart Hb (HB-201) & Test Strip

- Results in less than 10 sec.
- Results with HCT value.
- Low sample volume 10 µl.
- Measuring range 4.5 - 25.6 g/dL.
- Large memory of 1000 tests.
- Auto power off when not in use.
- Accurate and reliable result & Large display size.
- Data transmission to PC for Data management (optional).



Ordering Information:

Product Code	Product Name
SHB201QB	Smart Hb Meter (HB-201)
SHS050QB	Smart Hb Test Strip (HS-101)

Electrolyte Analyzers

InnoLyte

InnoLyte
Electrolyte Analyzer

- Based on ISE (Ion Selective Electrode Direct Measurement).
- 5.6-inch colour touch screen display.
- Throughput: 60 samples/hour.
- Sample type - Serum, others.
- Up-to five parameter options starting from K⁺ & Na⁺.
- Calibration: Automatic & manual.
- Closed reagent system with real time monitoring.
- In-built thermal printer.
- Up to 50,000 results storage.
- RS - 232 interface & LIS support.
- USFDA Listed & CE certified.

Ordering Information:

Product Code	Product Name
GE3001PS	InnoLyte (K ⁺ , Na ⁺ Parameters)

Innolyte Plus

- Based on ISE (Ion Selective Electrode Direct Measurement).
- 7.0 colour LCD touch screen.
- Throughput - 60 samples per/hour.
- Sample type - serum, plasma, whole blood and urine.
- pH corrected Ca^{++} Values.
- Automatic & manual calibration mode.
- Electrode shielding system.
- Alphanumeric sample ID input & connect bar-code scanner.
- Automated sampler (optional).
- RS - 232 interface & LIS support.
- More than 100,000 test results can be stored.
- CE certified.



InnoLyte^{Plus}
Electrolyte Analyzer

Ordering Information:

Product Code	Product Name
XI931AT	Innolyte Plus (K^+ , Na^+ Parameters)

Urine Chemistry Analyzers

Uriscan[®] Optima



URISCAN[®] Optima
World's Fastest Urine Chemistry Analyzer

- Based on reflectance photometer method.
- 4.3" Touch screen TFT LCD.
- Throughput 36 or 300 tests/ in two different modes.
- Workable with 2 - 11 parameter strips.
- Memory capacity of 2000 tests results.
- In-built thermal printer.
- 2 USB ports & RS-232 for interface.
- USFDA & CE certified.

Ordering Information:

Product Code	Product Name
M7YD	Uriscan Optima

Uriscan[®] Pro

- Based on reflectance photometer method.
- 7.0" Touch screen TFT LCD.
- Throughput 720 samples/ hour.
- Workable with 2-11 parameter strips.
- Memory capacity of 3000 tests results, 1000 controls & 30 calibrations.
- Three level QC software function.
- In-built thermal printer & external barcode reader.
- RS-232C & USB interface.
- USFDA & CE certified.



URISCAN[®] PRO
World's Fastest Urine Chemistry Analyzer

Ordering Information:

Product Code	Product Name
M9YD	Uriscan Pro

UriScan® Super⁺

- Random access fully automated urine chemistry analyzer.
- 8.4" large colour touch screen LCD.
- Throughput -200 tests/hour.
- Measurement method -CCD image scanner, refractometer & turbidimeter.
- 100 sample loading capacity with STAT facility.
- Microalbumin & creatinine in 13 parameter Super⁺ strip cassette.
- Unique cassette type reagent supply system.
- Integrated barcode reader.
- Interfacing port 4USB, RS-232C, serial interface, 2PS/2, parallel, ethernet ports for LIS.
- CE certified.




URISCAN® Super⁺
High-end Fully Automated Urine Chemistry Analyzer

Ordering Information:

Product Code	Product Name
M8YD	Uriscan Super Plus

ELISA Readers

Stat Fax® 4700



Stat Fax® 4700
ELISA Strip Reader

Ordering Information:

Product Code	Product Name
SF4700AT	Elisa Microstrip Reader

- Microwell strip reader based on bi-chromatic photometric measurements.
- Intuitive interface with touch screen display.
- Linear Measurement Range: 0.0 to 3.0 Abs Units.
- Three strip loading capacity.
- Reads 12 - well strip in 30 seconds.
- In-built printer.
- Standard four filter model (405 nm, 450 nm, 492 nm, 630 nm), other options available.
- User programmable test menu & various calculation modes.
- USFDA & CE certified

Stat Fax® 4200

- Compact standalone 8 – channel microplate reader.
- 5.7" touch-screen display.
- Measures absorbance of 96 wells in 12 seconds.
- Linear measurement range: 0.0 to 4.0 abs Unit.
- Four filter model (405 nm, 450 nm, 492 nm, 630 nm), other options available.
- Reading mode monochromatic & bi-chromatic.
- Reading mode monochromatic & bi-chromatic
- Open system with selectable plate formatting.
- In- Built printer.
- USFDA & CE certified



Stat Fax® 4200
ELISA Microplate Reader

Ordering Information:

Product Code	Product Name
SF4200AT	Elisa Microplate Reader

QL-9600



Q-lisa® Reader
QL-9600

- Based on bi - chromatic photometric measurement.
- 8-channel optical fibre system.
- Touch screen, large LCD display.
- Linear measurement range: 0.0 to 3.5 abs Units.
- Windows operation software.
- 5 second reading time for the plate.
- Four filter model (405 nm, 450 nm, 492 nm, 630 nm).
- Storage of up to 100 programs, 10,000 sample records.
- Multi-assay feature enables up to 12 different assays.
- Adjustable plate shaking speed & time.
- Powerful QC function, levey jennings plot graph.
- RS-232,USB & SD card interface.
- CE Certified.

Ordering Information:

Product Code	Product Name
QLR9600QB	Microplate Reader

ELISA Washers

Stat Fax® 2600

- 8-probe head micro plate washer.
- Two line LCD display with keyboard.
- Washes flats & round bottom wells.
- Automatic precision washing of microplates and microstrips .
- Full plate wash process time 55 seconds.
- Automatic calibration alignment and last row detection.
- Four filter model (405 nm, 450 nm, 492 nm, 630 nm).
- Six factory programmed wash/rinse modes.
- Memory 50 wash protocols.
- USFDA & CE certified.



Stat Fax® 2600
Microwell Plate Washer

Ordering Information:

Product Code	Product Name
SF2600AT	Microplate Washer

QL-4800



Q-lisa® Washer
QL-4800

- 8-way manifold, automatic and manual positioning.
- User Friendly operations system with large LCD display.
- Automatic monitoring for vacuum and pressure, automatic rinse cycle.
- Washes flat, U and V-bottom strips and plate.
- Large memory of up to 100 wash protocols.
- CE certified

Ordering Information:

Product Code	Product Name
QLW4800QB	Microplate Washer

Biomedical Systems

VAPRO® 5600

- Based on the vapor pressure depression principle.
- Capable of routine osmolality determinations on whole blood, biological samples, viscous liquids, thin slices of plant & animal tissue etc.
- Osmolality measured at room temperature without changing natural state.
- Self-cleaning thermocouple process reduces maintenance.
- Sample volumes from 2 to 10 µL.
- Test result in 90 seconds.
- Measuring range from 20 to 3200 mmol/kg.
- RS-232 (ASCII format) & USB – Slave.



Vapro

Ordering Information:

Product Code	Product Name
MODEL5600ET	Vapro Vapor Pressure Osmometer

Macroduct®


**MACRODUCT®
ADVANCED**
Sweat Induction and Collection

- Stimulation & collection of sweat for the diagnosis of cystic fibrosis.
- Easy to use touch screen interface.
- Flexible data management options.
- Meets protocols and traceability requirements.
- Supplied with Macroduct® Advanced Start-up Kit

Ordering Information:

Product Code	Product Name
MODEL3710ET	Macroduct Adv. Sweat Collection System

Sweat-Chek


- Sweat analysis for the diagnosis of cystic fibrosis.
- Designed specifically for use with the Macroduct® advanced sweat collector.
- It can measure the electrolyte concentration of a sweat specimen as small as 6 to 10 microliters.
- Reportable measurement range is 0 to 150 mmol/L.
- Supplied with Sweat-Chek™ Start-up Kit.


Sweat-Chek™
Sweat Conductivity Analyzer
Ordering Information:

Product Code	Product Name
MODEL3120ET	Sweat-Chek Sweat Conductivity Analyzer

Nanoduct®

- Complete integrated system for inducing and analysing sweat for cystic fibrosis (CF) diagnosis.
- Initial sweating rate displayed automatically.
- Requires only 3 µL of sweat.
- Shorter stimulation Time.
- Neonate-to-adult capability.
- Automatic averaging of conductivity.
- Compact Handy device.
- Supplied with Nanoduct® start-up kit and lab report software.



Ordering Information:

Product Code	Product Name
NANODUCTSYS1030ET	Nanoduct (Neonatal Sweat Analysis System)

ChloroChek®



**CHLOROCHek®
CHLORIDOMETER®**

- The Chloridometer® is a coulometric titrator designed to determine chloride ion concentrations in sweat samples .
 - Required 10 microliters of sweat sample.
 - Results in less than 20 seconds.
 - Intuitive touchscreen interface.
 - Results can be exported for complete traceability.
- Controls are available, with three levels, to ensure the highest quality control procedures to meet regulatory requirements.
- Supplied with ChloroChek® start-up Kit

Ordering Information:

Product Code	Product Name
MODEL3400ET	ChloroChek Chloridometer

Nucleic Acid Extraction Systems

Q-Xtra^{Plus}

- Fully automatic based on magnetic bead purification technology.
- Isolates pure, high quality DNA or RNA from variety of samples.
- Throughput up to 32 samples in nine minutes.
- Collection efficiency (Recovery Rate) $\geq 98\%$.
- Available with pre-filled reagent kits.
- Built – in UV sterilization facility.
- Five pre-set programmes and 500 programmable options.
- USFDA Listed & CE certified.



Q-Xtra^{Plus}

Ordering Information:

Product Code	Product Name
QXPM32PS	Q-Xtra Plus

Q-Xtra 96



Q-Xtra 96

- Automated system based on magnetic bead purification technology.
- Throughput of 96 samples processing in 15 minutes.
- Isolate DNA, RNA, proteins and cells from a variety of samples.
- Can store more than 100 user definable programs.
- Collection efficiency (Recovery Rate) $\geq 95\%$.
- Built-in UV light Sterilization facility.
- USFDA Listed & CE certified

Ordering Information:

Product Code	Product Name
QXPM96PS	Q-Xtra 96

REAGENTS AND CONSUMABLES

Q-line®



Q-line® S⁺

Q-line® PLATINUM



Q-lisa®

Q-line® Rapid



Q-line® Molecular



Q-Line® Generic Clinical Chemistry Liquid Stable Reagents



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
Albumin	ALBG0100QB	End Point - BCG	2 x 50 mL
	ALBG0250QB		2 x 125 mL
ALP (DEA)	ALPG0050QB	Kinetic, DGKC/SCE (DEA Buffer)	50 mL
	ALPG0100QB		2 x 50 mL
	ALPG0550QB		5 x 50 mL
ALP (IFCC)	ALPIG050QB	Kinetic, IFCC (AMP Buffer)	50 mL
	ALPIG100QB		2 x 50 mL
	ALPIG250QB		5 x 50 mL
ALT/GPT	ALTG0050QB	Kinetic UV IFCC w/o P5P	50 mL
	ALTG0100QB		2 x 50 mL
	ALTG0550QB		5 x 50 mL
Amylase	AMLG0050QB	Kinetic , Substrate CNPG3	2 x 25 mL
	AMLG0100QB		2 x 50 mL
AST/GOT	ASTG0050QB	Kinetic UV IFCC w/o P5P	50 mL
	ASTG0100QB		2 x 50 mL
	ASTG0550QB		5 x 50 mL
Bilirubin Direct	BIDG0250QB	End Point - Modified Malloy-Evelyn	2 x 125 mL
Bilirubin Total	BITG0250QB		2 x 125 mL
Bilirubin (T&D)	BTDG0250QB		2 x 125 mL
Bilirubin (T&D) Factor Based	BTDFS0250QB	End Point - Diazo Method	2 x 125 mL
Calcium	CALG0100QB	End Point - Arsenazo	2 x 50 mL
	CALG0250QB		2 x 125 mL
Cholesterol	CHOG0100QB	End Point - Enzymatic (CHOD/PAP)	2 x 50 mL
	CHOG0300QB		6 x 50 mL
	CHOG0600QB		6 x 100 mL
Chloride	CHLG0100QB	End Point - Thiocyanate	2 x 50 mL
Creatinine	CRTG0100QB	Jaffe Kinetic	2 x 50 mL
	CRTG0250QB		2 x 125 mL
Glucose PAP	GLUG1000QB	End Point - Enzymatic (GOD/ PAP)	2 x 500 mL
Gamma-GT	GGTG0050QB	IFCC Method	50 mL
HDL Direct	HDLG0080QB	Enzymatic - Direct Homogenous	80 mL with cal
	HDLG0380QB		380 mL
LDL Direct	LDLG0080QB	Enzymatic - Direct Homogenous	80 mL with cal
	LDLG0380QB		380 mL
HDL/LDL Calibrator	HDL/LDL- CAL	Human Serum Based	1 x 1 mL
LDH-P	LDPG0050QB	Kinetic. UV. - Pyruvate	50 mL
LIPASE	LPSG0020QB	Enzymatic Kinetic	20 mL with cal
	LPSG0040QB		2 x 20 mL with cal
	LPSGWC40QB		2 x 20 mL
Magnesium	MGXG0100QB	End Point - Xylidyl Blue	2 x 50 mL
Phosphorus	PHOG0100QB	End point UV - Phosphomolybdate	2 x 50 mL
Total Protein	TPRG0100QB	End Point - Biuret	2 x 50 mL
	TPRG0250QB		2 x 125 mL

Q-Line® Generic Clinical Chemistry Liquid Stable Reagents



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
Triglycerides	TRGG0100QB	End Point - Enzymatic (GPO/PAP)	2 x 50 mL
	TRGG0300QB		6 x 50 mL
Urea UV	UREG0050QB	Kinetic - Enzymatic UV (GLDH)	50 mL
	UREG0100QB		2 x 50 mL
	UREG0550QB		5 x 50 mL
Uric Acid	URIG0100QB	End Point - Enzymatic PAP (Uricase)	2 x 50 mL
	URIG0300QB		6 x 50 mL
MICROWASH (Innolab 200)	CSIL0300QB	Flowcell Cleaning Solution (Innolab 200)	6 x 50 mL
MICROWASH (Microlab 300)	CSML0300QB	Flowcell Cleaning Solution (Microlab 300)	6 x 50 mL

Q-Line® - Turbilatex Reagents



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
CRP Turbilatex	CRPG0050QB	Turbilatex	50 mL

LABKIT® Biochemistry Reagents from "Chemelex, Spain"



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
α - Amylase	30140LB	Kinetic-CNPG3	2 x 50 mL
Acid Phosphatase	30110LB	Kinetic - Colorimetric Hillmann	19 x 2 mL
Albumin	30120SLB	Bromocresol Green	2 x 50 mL
	30120LB		2 x 125 mL
ALP	30133SLB	PNPP (DEA Buffer)	50 mL
	30133S2LB		2 x 50 mL
	301333LB		125 mL (100 + 25)
Bilirubin Direct	30150LB	End Point - DMSO	2 x 125 mL
Bilirubin T + D	30157SLB	End Point - DMSO	2 x 50 mL
	30157LB		2 x 125 mL
Bilirubin Total	30155LB	End Point - DMSO	2 x 125 mL
Calcium	30165SLB	End Point - OCPC	2 x 50 mL
	30165LB		2 x 125 mL
Calcium Arsenazo III	30160SLB	End Point - Aarsenazo III	2 x 50 mL
	30160LB		2 x 125 mL
Chloride	30170LB	Thiocyanate	2 x125 mL
Cholesterol	30183SLB	End Point - CHOD-POD	2 x 50 mL
	30183LB		2 x125 mL
	30180LB		9 x 50 mL (Lyo)
Cholinesterase	30190LB	Kinetic - Butyrylthiocholine	20 x 3 mL
CK-MB LQ	30204SLB	Kinetic - UV Immunoinhibition	25 mL
CK-MB Control	30947LB	Human Serum Based	3 mL

LABKIT® Biochemistry Reagents from "Chemelex, Spain"



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
CK-NAC LQ	30201SLB	Kinetic - UV	25 mL
Copper	30205LB	Color 3-5 Di Br-PAESA	5 x 10 mL
Creatinine	30210/LB	Kinetic-Jaffe	2 x 50 mL
	30210LB		2 x 125 mL
Citric Acid	30169LB	Enzymatic	5 x 20 mL
Citric / Oxalic Ac. CTROL	30948LB	-	2 x 5 mL
GGT	30223SLB	Carboxy substrate	50 mL
GOT (AST) - LQ	30243SLB	Kinetic - Enzymatic	50 mL
	30243S2LB		2 x 50 mL
	30244LB		2 x (100 + 25) mL
GPT (ALT) - LQ	30253SLB	Kinetic - Enzymatic	50 mL
	30253S2LB		2 x 50 mL
	30254LB		2 x (100 + 25) mL
HDL Cholesterol Direct	30187LB	Direct - Enzymatic	80 mL with cal
Iron	30270LB	End Point - Ferrozine	4 x 50 mL
Labtrol H Calibrator	30960LB	Human Serum Based	10 x 3 mL
Labtrol H Normal	30950LB	Human Serum Based	4 x 5 mL
Labtrol H Pathological	30955LB	Human Serum Based	4 x 5 mL
Lactate	30280LB	Enzymatic LO - POD.	10 x 10 mL (Lyo)
LDH-P LQ	30293SLB	Kinetic Pyruvate (DGKC)	50 mL
LDL Cholesterol Direct	30185LB	Direct-Enzymatic	40 mL with cal
Lipase	30297LB	Enzymatic - Kinetic	4 x 10 mL
Magnesium	30300LB	End Point - Calmagite Method	2 x 125 mL
Oxalic Acid	30305LB	Enzymatic - Oxalate oxidase	20 Tests
Phospholipids	30320LB	Enzymatic - CHO-POD.	10 x 10 mL
Phosphorus UV	30330LB	End Point - UV Phosphomolybdate	2 x 125 mL
TIBC	30340LB	Saturation Precipitation	50 Tests
Total Proteins	30350SLB	End Point - Biuret	2 x 50 mL
	30350LB		2 x 125 mL
Total Proteins Urine & CSF (Microprotein)	30310LB	End Point - Pyrogallol Red	2 x 125 mL

LABKIT® Biochemistry Reagents from "Chemelex, Spain"


PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
Triglyceride	30364SLB	End Point - GPO-PAP	2 x 50 mL
Urea UV - LQ	30375SLB	Kinetic - Enzymatic UV-GLDH	50 mL
	30375S2LB		2 x 50 mL
	30375LB		2 x (100 + 25) mL
Uric Acid	30393SLB	End Point - Enzymatic-Uricase Method	2 x 50 mL
	30391LB		9 x 50 mL (Lyo)
Zinc	30400LB	5-Br-PAPS-Color	5 x 10 mL

LABKIT® Biochemistry Reagents Serology Range (Slide-Agglutination)


PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
ASO Latex Kit	40100LB	Slide Agglutination	50 Test.
	40101LB		100 Test.
CRP Latex Kit	40110LB	Slide Agglutination	50 Test.
	40111LB		100 Test.
RF Latex Kit	40110LB	Slide Agglutination	50 Test.
	40111LB		100 Test.

LABKIT® Biochemistry Turbilatex Reagents


PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
ASO Turbi with CAL	41015LB	Turbilatex	50 mL
RF Turbi with CAL	41025LB	Turbilatex	50 mL
CRP Turbi with CAL	41035LB	Turbilatex	50 mL
ASO-CRP-RF Control High	41012LB	Turbilatex	4 x 1 mL (Lyo)
ASO-CRP-RF Control Low	41022LB	Turbilatex	4 x 1 mL (Lyo)
Ferritin	41065LB	Turbilatex	50 mL with Cal
Ferritin Control	41062LB	Turbilatex (Liquid stable)	2 x 1 mL
Microalbumin	41095LB	Turbilatex	50 mL with Cal
Microalbumin control	41092LB	Turbilatex	1 x 2 mL
Lp (a) Turbilatex	41040LB	Turbilatex	52 mL with Cal
Lp (a) Control	41042LB	Turbilatex	1 x 1 mL
β2-Microglobulin	41055LB	Turbilatex	50 mL with Cal
β2-Microglobulin Control	41052LB	-	1 x 2 mL
Immunoglobulin IgE	41070LB	Turbilatex	18 (R1=13 + R2=5) mL
IgE Calibrator Set	41071LB	Turbilatex (Liquid stable)	6 x 1 mL
IgE Control Set	41072LB	Turbilatex (Liquid stable)	2 x 1 mL
α1- Antitrypsin	42031LB	Immunoturbidimetry	50 mL

LABKIT Biochemistry Immunoturbidimetry Reagent



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
α 1-acid Glycoprotein	42051LB	Immunoturbidimetry	50 mL
Apolipoprotein A1	42011LB	Immunoturbidimetry	50 mL
Apolipoprotein B	42021LB	Immunoturbidimetry	50 mL
Calibrator DUAL APO A1/B	42230LB	Immunoturbidimetry	1 x 1 mL
Control Serum DUAL APO A1/B	42240LB	Immunoturbidimetry	1 x 1 mL
Complement C3	42091LB	Immunoturbidimetry	50 mL
Complement C4	42101LB	Immunoturbidimetry	50 mL
Antithrombin III	42061LB	Immunoturbidimetry	50 mL
Ceruloplasmin	42071LB	Immunoturbidimetry	50 mL
Haptoglobin	42111LB	Immunoturbidimetry	50 mL
Immunoglobulin IgA	42121LB	Immunoturbidimetry	50 mL
Immunoglobulin IgG	42131LB	Immunoturbidimetry	50 mL
Immunoglobulin IgM	42141LB	Immunoturbidimetry	50 mL
Pre-Albumin	42191LB	Immunoturbidimetry	50 mL
Transferrin	42201LB	Immunoturbidimetry	50 mL
General Calibrator Immunoproteins	42210LB	"Immunoturbidimetry Assay Calibrator for 13 parameters (Liquid stable)"	1 x 2 mL
Control Serum Immunoproteins	42220LB	"Immunoturbidimetry Assay Control for 13 parameters (Liquid stable)"	4 x 1 mL

Q-Line® S⁺ Selectra System Pack Biochemistry Reagents

PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
ALBUMIN SP	ALBS0240QB	End Point - BCG	12 x 20 mL
ALP (DEA) SP	ALPS0100QB	Kinetic - DGKC/SCE (DEA Buffer)	4 x 25 mL
ALP (IFCC) SP	ALPIS200QB	Kinetic - IFCC (AMP Buffer)	8 x 25 mL
ALT/GPT SP	ALTS0200QB	Kinetic UV - IFCC without P5P	8 x 25 mL
AMYLASE SP	AMLS0120QB	Kinetic - Substrate CNPG3	6 x 20 mL
AST/GOT SP	ASTS0200QB	Kinetic UV - IFCC without P5P	8 x 25 mL
BILIRUBIN DIRECT SP	BIDS0200QB	End Point - Modified Malloy-Evelyn	8 x 25 mL
BILIRUBIN TOTAL SP	BITS0200QB	End Point - Modified Malloy-Evelyn	8 x 25 mL
CALCIUM ARSENAZO SP	CALS0240QB	End Point - Arsenazo III	12 x 20 mL
CHLORIDE SP	CHLS0240QB	End Point - Thiocyanate	12 x 20 mL
CHOLESTEROL SP	CHOS0240QB	End Point - Enzymatic (CHOD/PAP)	12 x 20 mL
CREATININE PAP SP	CRES0224QB	Kinetic - Enzymatic - PAP	8 x 28 mL

Q-Line® S+ Selectra System Pack Biochemistry Reagents



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
GAMMA-GT SP	GGTS0200QB	Kinetic - IFCC	8 x 25 mL
GLUCOSE PAP SP	GLUS0240QB	End Point - Enzymatic (GOD/ PAP)	12 x 20 mL
LDH-P SP	LDPS0100QB	Kinetic UV - IFCC	4 x 25 mL
LIPASE SP	LPSS0100QB	Kinetic - Enzymatic	4 x 25 mL
MAGNESIUM SP	MGXS0096QB	End Point - Xylidyl Blue	12 x 8 mL
PHOSPHORUS SP	PHOS0120QB	End point - UV Phosphomolybdate	6 x 20 mL
TOTAL PROTEIN SP	TPRS0240QB	End Point - Biuret	12 x 20 mL
TRIGLYCERIDE SP	TRGS0240QB	End Point - Enzymatic (GPO/PAP)	12 x 20 mL
UREA UV SP	URES0200QB	Kinetic - Enzymatic UV (GLDH)	8 x 25 mL
URIC ACID SP	URIS0240QB	End Point - Enzymatic PAP (Uricase)	12 x 20 mL
ELICAL 2	CALI0550ET	Multiparametric Calibrator (Human Serum Based)	4 x 3 mL
ELITROL I	CONT0060ET	Normal Multiparametric Control (Human Serum Based)	10 x 5 mL
ELITROL II	CONT0160ET	Abnormal Multiparametric Control (Human Serum Based)	10 x 5 mL
System Solution	SLSY5905ET	-	1000 mL
	SYSL1000QB		1000 mL

ASTIRA Series Chemistry System Packs Biochemistry Reagents



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
ALBUMIN ASP (C-MS)	ASALBU0500QB	End Point Bromocresol Green (BCG)	500 Tests
ALBUMIN ASP (C-MS/CX400)	ASALBU1000QB	End Point Bromocresol Green (BCG)	1000 Tests
ADA Kit ASP (C-MS/CX400)	ASADAE0200QB	Enzymatic	200 Tests
ALP (DEA) ASP (C-MS)	ASALPD0330QB	Kinetic Colorimetric DGKC (DEA Buffer)	330 Tests
ALP (DEA) ASP (C-MS/CX400)	ASALPD0660QB	Kinetic Colorimetric DGKC (DEA Buffer)	660 Tests
ALP (AMP) ASP (C-MS)	ASALPI0330QB	Kinetic Colorimetric IFCC (AMP Buffer)	330 Tests
ALP (AMP) ASP (C-MS/CX400)	ASALPI0660QB	Kinetic Colorimetric IFCC (AMP Buffer)	660 Tests
ALT/GPT ASP (C-MS)	ASALAT0530QB	Kinetic UV IFCC without P5P	530 Tests
ALT/GPT ASP (C-MS/CX400)	ASALAT1000QB	Kinetic UV IFCC without P5P	1000 Tests
AST/GOT ASP (C-MS)	ASASAT0530QB	Kinetic UV IFCC without P5P	530 Tests
AST/GOT ASP (C-MS/CX400)	ASASAT1000QB	Kinetic UV IFCC without P5P	1000 Tests
AMYLASE ASP (C-MS/CX400)	ASAMYL0200QB	Kinetic Colorimetric Substrate CNPG3	200 Tests
BILIRUBIN DIRECT ASP (C-MS)	ASBIDI0530QB	End Point Modified Malloy-Evelyn	530 Tests
BILIRUBIN DIRECT ASP (C-MS/CX400)	ASBIDI1000QB	End Point Modified Malloy-Evelyn	1000 Tests
BILIRUBIN TOTAL ASP (C-MS)	ASBITO0530QB	End Point Modified Malloy-Evelyn	530 Tests
BILIRUBIN TOTAL ASP (C-MS/CX400)	ASBITO1000QB	End Point Modified Malloy-Evelyn	1000 Tests
BILIRUBIN TOTAL ASP (C-MS/CX400)	ASBITO1000QB	End Point Modified Malloy-Evelyn	1000 Tests
CALCIUM ARS III ASP (C-MS/CX400)	ASCALA0200QB	End Point (Arsenazo III)	200 Tests
CHLORIDE ASP (C-MS/CX400)	ASCHLO0150QB	End Point Colored Complex (Thiocyanate)	150 Tests
CHOLESTEROL ASP (C-MS/CX400)	ASCHOL0500QB	End Point Enzymatic (CHOD/PAP)	500 Tests
CREATININE PAP ASP (C-MS)	ASCRTE0500QB	Kinetic Enzymatic - PAP	500 Tests
CREATININE PAP ASP (C-MS/CX400)	ASCRTE1000QB	Kinetic Enzymatic - PAP	1000 Tests

ASTIRA Series Chemistry System Packs Biochemistry Reagents



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
CRP Turbilatex ASP (C-MS/CX400)	ASCRPQ0250QB	Turbilatex	250 Tests
GAMMA-GT ASP (C-MS/CX400)	ASGGTI0330QB	Kinetic - IFCC	330 Tests
G-6-PD ASP (C-MS/CX400)	ASG6PD0150QB	UV Kinetic	150 Tests
GLUCOSE PAP ASP (C-MS/CX400)	ASGLUC1000QB	End Point Enzymatic (GOD/ PAP)	1000 Tests
HbA1c Direct ASP (C-MS/CX400)	ASHBAC0140QB	Latex Immunoturbidimetry	140 Tests
HDL Cholesterol Direct ASP (C-MS)	ASHDLC0400QB	Direct Homogenous Enzymatic	400 Tests
HDL Cholesterol Direct ASP (C-MS/CX400)	ASHDLC0800QB	Direct Homogenous Enzymatic	800 Tests
LDL Cholesterol Direct ASP (C-MS)	ASLDLC0400QB	Direct Homogenous Enzymatic	400 Tests
LDL Cholesterol Direct ASP (C-MS/CX400)	ASLDLC0800QB	Direct Homogenous Enzymatic	800 Tests
LDH-P ASP (C-MS/CX400)	ASLDHP0330QB	Kinetic UV IFCC	330 Tests
LIPASE ASP (C-MS/CX400)	ASLIPS0160QB	Enzymatic Kinetic	160 Tests
MAGNESIUM (XYLIDYL) ASP (C-MS/CX400)	ASMGXB0150QB	End Point - Xylidyl Blue	150 Tests
PHOSPHORUS ASP (C-MS/CX400)	ASPHOS0150QB	Phosphomolybdate UV End Point	150 Tests
TOTAL PROTEIN ASP (C-MS/CX400)	ASTPRT0450QB	End Point Biuret	450 Tests
TRIGLYCERIDES ASP (C-MS/CX400)	ASTRIG0500QB	End Point Enzymatic (GPO/PAP)	500 Tests
UREA/BUN UV ASP (C-MS)	ASUREA0400QB	Kinetic Enzymatic UV (GLDH)	400 Tests
UREA/BUN UV ASP (C-MS/CX400)	ASUREA0750QB	Kinetic Enzymatic UV (GLDH)	750 Tests
URIC ACID ASP (C-MS)	ASURIC0330QB	End Point Enzymatic PAP (Uricase)	330 Tests
URIC ACID ASP (C-MS/CX400)	ASURIC1060QB	End Point Enzymatic PAP (Uricase)	1060 Tests
Anti-Bacterial Phosphor Free Detergent (C-MS)	ASABPD0040QB	-	2 x 20 mL
Anti-Bacterial Phosphor Free Detergent (C-MS/CX400)	ASABPD0500QB	-	500 mL
Alkaline Detergent (C-MS/CX400)	ASALKD1000QB	-	1000 mL
Alkaline Detergent (C-MS/CX400)	ASALKD2000QB	-	2000 mL

Q-Line® PLATINUM Special Biochemistry Reagents



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
ADA Kit	ADAG0030QB	Enzymatic	30 mL
	ADAG0075QB		75 mL
ADA Calibrator	ADACAL01QB	Lyophilized	1 x 1 mL
ADA Control	ADACON02QB	Lyophilized	2 x 1 mL
ACE Kit	ACES0020QB	FAPGG Substrate	2 x 10 mL
ACE Calibrator	ACECAL01QB	Lyophilized	1 x 1 mL
ACE Control	ACECON02QB	Lyophilized	2 x 1 mL
D-Dimer Kit	DDIS0028QB	Latex Immunoturbidimetry	1 x 28 mL
D-Dimer Calibrator	DDICAL05QB	Lyophilized	5 x 1 mL
D-Dimer Controls	DDICON02QB	Lyophilized	2 x 1 mL
HbA1c Direct Kit	HBAS0024QB	Latex Immunoturbidimetry	1 x 24 mL
HbA1c Direct Calibrator	HBACAL05QB	Lyophilized	5 x 1 mL
HbA1c Direct Control	HBACON02QB	Lyophilized	2 x 1 mL

Q-Line® PLATINUM Special Biochemistry Reagents



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
Homocysteine Kit	HCYS0025QB	Enzyme Cycling	1 x 25 mL
Homocysteine Calibrator	HCYCAL05QB	Liquid stable	5 x 1 mL
Homocysteine Control	HCYCON02QB	Liquid stable	2 x 1 mL
G-6-PD	G6PDS024QB	UV Kinetic	1 x 24 mL
G-6-PD Control	G6PCON04QB	----	2 x 2 mL
IgE Kit	IGEG0032QB IGEG0064QB	Latex Immunoturbidimetry	1 x 32 mL 1 x 64 mL
IgE Controls	IGECON02QB	Latex Immunoturbidimetry	2 x 1 mL
Vitamin- D Kit	VTDS0025QB	Latex Immunoturbidimetry	1 x 25 mL
Vitamin- D Calibrator	VTDCAL05QB	(Liquid)	5 x 1 mL
Vitamin-D Controls	VTDCON02QB	(Liquid)	2 x 1 mL
PCT Kit	PCTS0024QB	Latex Immunoturbidimetry	1 x 24 mL
PCT Calibrator	PCTCAL06QB	Six-Point Calibration (Liquid)	6 x 1 mL
PCT Control	PCTCON02QB	(Liquid)	2 x 1 mL

Clinical Chemistry Reagents From DIAZYME, USA



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
Adenosine Deaminase kit	DZ117A-KY1	Enzymatic (Liquid Stable Dual Reagent)	75 mL (R1:1x50 mL+R2:1x25mL)
ADA Calibrator	DZ117A-CAL	Lyophilized	1x1 mL
ADA Control	DZ117A-CON	Lyophilized	2x1 mL
Direct Enzymatic HbA1c Kit	DZ168A-K	Enzymatic Single Channel (Liquid Stable-3 Reagensts)	34 mL (R1a-16.8 mL, R1b-7.2 mL, R2-10 mL Lysis Buffer-1x30 mL)
HbA1c Calibrator	DZ168A-CAL	Lyophilized	2 x 0.5 mL
HbA1c Control	DZ168A-CON	Lyophilized	2 x 0.5 mL
HDL Cholesterol Direct	DZ129A-KY1	PVS\PEGME Coupled	380 mL (R1- 5x56 mL, R2- 2x50 mL)
LDL Cholesterol Direct	DZ128A-KY1	PVS\PEGME Coupled	380 mL (R1- 5x56 mL, R2- 2x50 mL)
HDL/LDL Calibrator	DZ130A-CAL	Lyophilized	1 x 1 mL
HDL/LDL Control	DZ138A-CON	Lyophilized	3 x 3 mL
Total Bile Acids Assay	DZ042A-K DZ042A-K01	Enzyme Cycling	160 mL (R1:2 x 60 mL + R2:2 x 20 mL) 80 mL (R1:2 x 30 mL + R2:2 x 10 mL)
Total Bile Acids Assay Calibrator	DZ042A-CAL	Liquid Stable	1 x 2 mL
Total Bile Acids Control Set	DZ042A-CON	Lyophilized	2 x 1 mL
Cystatin C Assay	DZ133C-KY1	Immunoturbidimetric	72 mL (R1: 1 x 60 mL R2: 1 x 12 mL)
Cystatin C Calibrator Set	DZ133C-CAL	Liquid Stable	5 x 1 mL
Cystatin C Control Set	DZ133C-CON	Liquid Stable	2 x 1 mL
Human Kappa Free Light Chain Assay with Calibrator	DZ169A-KY1	Liquid Stable	29 mL (R1:1 x 22 mL+ R2:1 x 7 mL) + 5 x 1.5 mL (CAL)

Clinical Chemistry Reagents From DIAZYME, USA



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
Human Lambda Free Light Chain Assay with Cali	DZ170A-KY1	Liquid Stable	29 mL (R1:1 x 22 mL+ R2:1 x 7 mL) + 5 x 1.5 mL (CAL)
Kappa/Lambda Free Light Chain Control	DZ169A-CON	Liquid Stable	2 x 3 mL

Hemostasis Reagents From Stago & TCOAG, Ireland (A Stago Group Company)



PRODUCT NAME	PRODUCT CODE	PACK SIZE
Neoplastine® CI Plus 2 (ISI≈1.3)(PT)	00374ST	6 x 2 mL
Neoplastine® CI Plus 5 (ISI≈1.3) (PT)	00375ST	6 x 5 mL
CK PREST 2 (APTT)(Kaolin Activator)(Lyo)	00598ST	6 x 2 mL
STA CaCl ₂ (0.025 M)	00367ST	24 x 15 mL
STA Coag Control N+P (Lyo) (PT, APTT,FIB,TT,LUP)	00679ST	12 x 2 x 1 mL
Macro cuvettes for KC1 Delta™ with balls packed separately without ball dispenser	Z05100TC	1000 Nos.

Reagents & Consumables for Electrolyte Analyzer "INNOLYTE"



PRODUCT NAME	PRODUCT CODE	PACK SIZE
AB Standard Reagent (DS-I)(A:390+B:160)	IEAB0550PS	550 mL
Cleaning Solution	IECL0250PS	250 mL
Deprotein Solution	IEDS0020PS	20 mL
	IEDS0050PS	50 mL
	IEDS0125PS	125 mL
QC Solution (Low, Normal, High)	IECS0030PS	3 x 10 mL
Ref. Electrode Filler	IERF0010PS	10 mL
K ⁺ , Na, Cl ⁻ , Ca, Li Electrode Filler	IEEF0010PS	10 mL
PVC Electrode Activator Solution	902000103003402	250 mL
Pump Tube	P121532202200	1 Pcs.
Reagent Tube	P010910006300	1 Mtre

Reagents & Consumables for Electrolyte Analyser "INNOLYTE PLUS"



PRODUCT NAME	PRODUCT CODE	PACK SIZE
AB Standard Reagent	XI931FT Plus	1050 mL
Weekly Cleaning Solution	XI931-16-I-Plus	100 mL
Daily Cleaning Solution	XI931-15	100 mL
QC Solution (ampules)	XI931-14	50 x 2 mL
Urine Solution(Urine Diluent)	XI931-22-I-Plus	100 mL

Urine Chemistry Strips & Controls From Y.d. Diagnostics, Korea

Your Dream
Diagnose Everything Alive!

PRODUCT NAME	PRODUCT CODE	DESCRIPTION	PACK SIZE
URISCAN 2 GP strip	U20YD	Test parameters: Protein, Glucose	100 Strips
URISCAN 6 L strip	U25YD	Test parameters: Blood,Protein,Nitrite, Glucose, pH, Leucocytes	100 Strips
URISCAN 10 SGL strip	U39YD	Test parameters : Blood, Bilirubin, Urobilinogen, Ketones (Acetoacetic Acid), Protein,Nitrite, Glucose, pH, S.G (Specific Gravity), Leucocytes	100 Strips
URISCAN 11 strip	U41YD	Test parameters: Blood, Bilirubin, Urobilinogen, Ketones (Acetoacetic Acid), Protein,Nitrite, Glucose, pH,S.G, Leucocytes and Ascorbic Acid	100 Strips
ACR 02 Strip	U64YD	Albumin Creatinine & Albumin Creatinine Ratio	50 Strips
ACR 11 Strip	U63YD	Test parameters : Blood, Ketones (Acetoacetic Acid), Protein,Nitrite, Glucose, pH, S.G (Specific Gravity), Leucocytes, Ascorbic Acid, Microalbumin,Creatinine & ACR	100 Strips
URISCAN Super ⁺ Cassette	U81YD	Test parameters : Blood,Bilirubin, Creatinine, Urobilinogen, Ketones, Protein, Nitrite, Microalbumin, Glucose, pH, Leucocytes, Ascorbic acid, Color S.G (By Refractometer), Clarity(By Turbidimeter)	400 Strips
Uritrol - Urine Control (Level I,II,III)	U504YD	Assayed quality control urine to monitor and maintain the accuracy and precision of the urine strips and the analyzers	3 X10 ML
Rinse Solution	RINSYD	Probe Cleaning Solution for Uriscan Super	1 Litre

Q-Line[®] RAPID TEST Product Range

Q-line[®] Rapid

PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
HCG Test Card	HCGR0050QB	Immunochromatography	50 Test
Dengue Ab Test Card	DABR0010QB	Immunochromatography	10 Test
Dengue Antigen NS1 Test Card	DNSR0010QB	Immunochromatography	10 Test
Dengue Antigen NS1 Test Card (WB)	DNSWR010QB	Immunochromatography	10 Test
Dengue Combo(Ag & Ab) Test Card	DCMR0010QB	Immunochromatography	10 Test
Dengue Combo(Ag & Ab) Test Card	DCMR0025QB	Immunochromatography	25 test
Malaria Antigen (Pf/Pv) Test Card	MFVR0025QB	Immunochromatography	25 test
Malaria Antigen (Pf/Pv) Test Card	MFVR0050QB	Immunochromatography	50 Test
Malaria Antigen (Pf/Pv) P.O.C.Test Card	MFVPR001QB	Immunochromatography	1 Test
Malaria Antigen (Pf/Pan) Test Card	MFPR0025QB	Immunochromatography	25 tests
Malaria Antigen (Pf/Pan) Test Card	MFPR0050QB	Immunochromatography	50 Test
HbsAg Antigen Test Card (whole Blood/serum/plasma)	HBSR0025QB	Immunochromatography	25 test
HbsAg Antigen Test Card (whole Blood/serum/plasma)	HBSR0050QB	Immunochromatography	50 Test
HbsAg Antigen Test Card (serum/plasma)	HBSSR025QB	Immunochromatography	25 test
HbsAg Antigen Test Card (serum/plasma)	HBSSR050QB	Immunochromatography	50 Test
HIV 1/2 Ab Test Card (3rd Generation)	HIVR0025QB	Immunochromatography	25 Test
HIV 1/2 Ab Test Card (3rd Generation)	HIVR0050QB	Immunochromatography	50 Test
HIV-Ab/Ag (4th Gen) Test Card	HIVFR025QB	Immunochromatography	25 Test
HIV-Ab/Ag (4th Gen) Test Card	HIVFR050QB	Immunochromatography	50 Test
HCV Ab Test Card	HCVR0025QB	Immunochromatography	25 Test
HCV Ab Test Card	HCVR0050QB	Immunochromatography	50 Test
Syphilis Ab Card	SYLCR025QB	Immunochromatography	25 Test
HIV/Syphilis DUO Test Card	HIVSYR25QB	Immunochromatography	25 Test
HIV/Syphilis DUO Test Card	HIVSYR50QB	Immunochromatography	50 Test

Q-Lisa® ELISA Product Range



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
HIV -1/2 ELISA (4 th Genration)	HIVEF096QB	Microwell ELISA	96 Test
HIV -1/2 ELISA (3 rd Genration)	HIVE0096QB	Microwell ELISA	96 Test
HBsAg ELISA (3 rd Genration)	HBSE0096QB	Microwell ELISA	96 Test
HBsAg ELISA High Sensitive	HBSEH096QB	Microwell ELISA	96 Test
HCV ELISA (3 rd Genration)	HCVE0096QB	Microwell ELISA	96 Test
Covid-19 ELISA (IgG Antibody)	COVEIGG96QB	Microwell ELISA	96 Test
Dengue NS1 ELISA	DNSE0096B	Microwell ELISA	48 Test
Syphilis Antibody ELISA	SYPE0096QB	Microwell ELISA	96 Test

Q-Line® Molecular Diagnostics (MDX) -COVID Products



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
Viral Nucleic Acid Extraction	MDVNXS0050QB	Spin Column Method (Manual Application)	50 Test
Viral Nucleic Acid Extraction	MDVNXS0100QB	Spin Column Method (Manual Application)	100 Test
Viral DNA Extraction	MDVDXS0050QB	Spin Column Method (Manual Application)	50 Test
Viral DNA Extraction	MDVDXS0100QB	Spin Column Method (Manual Application)	100 Test
Viral RNA Extraction	MDVRXS0050QB	Spin Column Method (Manual Application)	50 Test
Viral RNA Extraction	MDVRXS0100QB	Spin Column Method (Manual Application)	100 Test
Bacterial Genomic DNA Extraction	MDBDXS0050QB	Spin Column Method (Manual Application)	50 Test
Bacterial Genomic DNA Extraction	MDBDXS0100QB	Spin Column Method (Manual Application)	100 Test
Viral Nucleic Acid Extraction	MDVNXM0050QB	Magnetic Bead Manual Method (Manual Application)	50 Test
Viral Nucleic Acid Extraction	MDVNXM0100QB	Magnetic Bead Manual Method (Manual Application)	100 Test
Viral Nucleic Acid Extraction	MDVN32Q032QB	Magnetic Bead Manual Method (Automated Application Q-xtra 32 Sample Extraction System)	32 Test
Viral Nucleic Acid Extraction	MDVN32Q096QB	Magnetic Bead Manual Method (Automated Application Q-xtra 32 Sample Extraction System)	96 Test
Viral DNA Extraction	MDVD32Q032QB	Magnetic Bead Manual Method (Automated Application Q-xtra 32 Sample Extraction System)	32 Test
Viral DNA Extraction	MDVD32Q096QB	Magnetic Bead Manual Method (Automated Application Q-xtra 32 Sample Extraction System)	96 Test
Viral RNA Extraction	MDVR32Q032QB	Magnetic Bead Manual Method (Automated Application Q-xtra 32 Sample Extraction System)	32 Test
Viral RNA Extraction	MDVR32Q096QB	Magnetic Bead Manual Method (Automated Application Q-xtra 32 Sample Extraction System)	96 Test
Bacterial Genomic DNA Extraction	MDBD32Q032QB	Magnetic Bead Manual Method (Automated Application Q-xtra 32 Sample Extraction System)	32 Test
Bacterial Genomic DNA Extraction	MDBD32Q096QB	Magnetic Bead Manual Method (Automated Application Q-xtra 32 Sample Extraction System)	96 Test
HBV RT-PCR Kit (Quantitative)	MDHBVT0025QB	RT-PCR	25 Test
HBV RT-PCR Kit (Quantitative)	MDHBVT0100QB	RT-PCR	100 Test
HCV RT-PCR Kit (Quantitative)	MDHCVT0025QB	RT-PCR	25 Test
HCV RT-PCR Kit (Quantitative)	MDHCVT0100QB	RT-PCR	100 Test
HIV RT-PCR Kit (Quantitative)	MDHIV1T025QB	RT-PCR	25 Test
HIV RT-PCR Kit (Quantitative)	MDHIV1T100QB	RT-PCR	100 Test
MTB RT-PCR Kit	MDMTBX0025QB	RT-PCR	25 Test
MTB RT-PCR Kit	MDMTBX0100QB	RT-PCR	100 Test
HBV/HCV/HIV RT-PCR Kit	MDHHHML025QB	RT-PCR	25 Test
HBV/HCV/HIV RT-PCR Kit	MDHHHML100QB	RT-PCR	100 Test
Dengue RT-PCR Kit	MDDNGU0025QB	RT-PCR	25 Test
Dengue RT-PCR Kit	MDDNGU0100QB	RT-PCR	100 Test

Hematology Reagents, Calibrators, Controls & ESR Tubes



PRODUCT NAME	PRODUCT CODE	PACK SIZE
Medonic M-Series Diluent (RFID)	1504460BM/IN1504460BM	20 Litre
Medonic M-Series Lyse (RFID)	1504461BM/IN1504461BM	5 Litre
Boule Control Low (M20/M32)	1504289BM	3 mL
Boule Control Normal (M20/M32)	1504288BM	3 mL
Boule Control High (M20/M32)	1504290BM	3 mL
Control (Low, Normal, High)(M20/M32)	1504022BM	6 x 4.5 mL
Human Calibrator (M20/M32)	1504025BM	1 x 3 mL
Medonic M51 Diluent	1504510BM	20 Litre
Medonic M51-L1 Lyse	1504511BM	200 mL
Medonic M51-L2 Lyse	1504512BM	500 mL
Boule Easy Cleaner	1504513BM	50 mL
Boule Cal-5 Diff A1(M51)	1504517BM	1 x 3 mL
Boule Con-5 Diff A1 Low (M51)	1504521BM	1 x 3 mL
Boule Con-5 Diff A1 Norm (M51)	1504520BM	1 x 3 mL
Boule Con-5 Diff A1 High (M51)	1504522BM	1 x 3 mL
Hematology Cleaning Kit (M20/M32)	CKMS1350PS	3 x 450 mL
Boule Cleaning Kit	1504111BM	3 x 450 mL
Sediscan ESR Vacuum Tube for Microsed/MIX RATE X20	PRD-PRV11B	1 Tube
Qualscan ESR Vacuum Tubes for Qualcyte 10 (With Smart Card)	ESRQS0200LBX	2 x 100 Tube

P.O.C. Products



PRODUCT NAME	PRODUCT CODE	METHOD	PACK SIZE
Smart Hb (HB-201)	SHB201QB	Digital Hemoglobin Meter	1 Pc
Smart Hb Test Strip (HS-101)	SHS050QB	Hb Strips for Smart Hb (HB-201)	1 Bottle (50 Strips)

Elitech Biomedical Systems Reagents



PRODUCT NAME	PRODUCT CODE	PACK SIZE
OPTIMOLE™ Osmolality Standard	OA-010ET	60 x 0.4 mL
OPTIMOLE™ Osmolality Standard	OA-029ET	60 x 0.4 mL
OPTIMOLE™ Osmolality Standard	OA-100ET	60 x 0.4 mL
OSMOCOLL® HNL Osmolality Controls	SS-0273ET	6 x 1mL
Sample Discs	SS-033ET	5000 nos
VAPRO® Desiccant Filter Cartridge	SS-238ET	1 nos.
Thermocouple Cleaning Solution	SS-003ET	60 mL
Thermocouple Cleaning Solution w/Anionic Detergent	SS-223ET	60 mL



Certificate of Registration

This certificate has been awarded to

Q-LINE BIOTECH PRIVATE LIMITED

Office—C-108, Naraina Industrial Area, Phase-1, Naraina, New Delhi, 110028, India

In recognition of the organization's Quality Management System which complies with

ISO 13485:2016

The scope of activities covered by this certificate is defined below

Design, Manufacture and Distribution of In-Vitro Diagnostic Reagents and Instruments

Certificate Number **84700/A/0003/UK/En**


A certificate number of 0001, confirms the Client has a single site Certified & the site is their Head Office or Main site in relation to the Certified scope with URS. A certificate number of 0002, or greater (e.g.: xxxx/M/0002/UK/En) refers to a client that has more than one site certified with URS, as such, the following statement shall apply - 'The validity of this certificate depends on the validity of the main certificate'.

Date of Issue of Certification Cycle	Issue Number	Certificate Expiry Date	Certification Cycle
19 October 2024	3	18 October 2027	3
Revision Date	Revision Number	Original Certificate Issue Date	Scheme Number
15 September 2024	0	19 October 2018	n/a

For detailed explanation for the data fields above, refer to <http://www.urs-holdings.com/logos-and-regulations>

Issued by

Mukesh Singhal - On behalf of the Schemes Manager



If there is any doubt as to the authenticity of this certificate, please do not hesitate to contact the Head Office of the Group on info@urs-certification.com. URS is a member of United Registrar of Systems (Holdings) Ltd, United House, 28 Poole Hill, Bournemouth, BH2 9PS, UK. Company Registration no. 5298466

1 / 1

Corporate Office :

Q-LINE BIOTECH LIMITED

298-281, Transport Nagar, Kanpur Road,
Adjacent Transport Nagar Metro Station,
Lucknow - 226012 (U.P), India
Phone - 0522-2433023


Branch Office :

KOLKATA

Fb161, Rajdanga Main Road, Kolkata ,
(West Bengal) - 700 107

ODISHA

Anand Bhawan, Plot No- 441/1621/1807, Divyanagar
Nuagaon, Side Lane Of Indian Oil Petrol Pump Nuagaon,
Bhubaneshwar, Odisha Khordha, Odisha - 751002



Q-LINE BIOTECH LIMITED

C -108, Naraina Industrial Area, Phase-1, New Delhi - 110028.

Ph: +91-11-45577407